

	CORE PELLETS					COATED PELLETS
	DRUG LAYERED CORE PELLETS		MATRIX CORE PELLETS			
structure of core pellets / coated pellets						
Glatt Technology	Wurster (bottomspray) 	CPS/rotor powderlayering 	CPS/rotor pelletisation 	Extrusion/Spheronisation 	MicroPx 	ProCell
process mode	batch	batch	batch	batch/continuous	continuous	
solid starting material	starter pellets ≥ 100 µm	starter pellets ≥ 100 µm	microcrystalline cellulose powder API etc.	microcrystalline cellulose powder API etc.	no	
processing of API/coating liquid solvent(s)	from a liquid (solution, suspension, emulsion, melt) water, organic solvent(s)	as powder water, organic solvent(s)	as powder water, organic solvent(s)	wetted mass of API + microcrystalline cellulose water, organic solvent(s)	from a liquid (solution, suspension, emulsion, melt) water, organic solvent(s)	
API dosage	strength / drug load		strength / drug load		strength / drug load	strength / drug load
low	••••• ~0,01 - 10%		na ~0,01 - 10%		••• ~0,01 - 10%	••• ~0,01 - 10%
medium	••••• ~10 - 50%		••• ~10 - 50%		••••• ~10 - 50%	••••• ~10 - 50%
high	••••• ~50 - 80%		••••• ~50 - 80%		••••• ~50 - 80%	••••• ~50 - 80%
very high	na ~80 - 100%		na ~80 - 100%		••••• ~80 - 100%	••••• ~80 - 100%
properties of API	chemically stable in water, organic solvents		sensitive to moisture		chemically stable in water, organic solvents, melt	
smallest pellet size achievable	~150 - 400 µm micropellets possible		~300 - 400 µm micropellets possible		~150 - 400 µm micropellets possible	
comment	preferred technology for drug layering with chemically stable APIs		preferred technology for medium / high dosed moisture sensitive APIs		high impact of physico-chemical properties of API on the process quality (sticking of wetted API)	often imperfect sphericity and surface smoothness potential impact on coating quality and performance
high potent / containment feasibility	smart high potent processing approach					usual drug load: 90 - 95% preferred technology for high drug loaded micropellets
						usual drug load: 90 - 95% preferred technology for high drug loaded pellets including temperature sensitive API
						smart high potent processing approach
						coating liquid: · solution · dispersion · melt
						solvents: · water · organic solvent(s)
						~200 - 500 µm coated micropellets possible
						preferred technology for all coating applications · controlled drug release · taste masking · seal coating · pH dependent coating · immediate release top coat
						smart high potent processing approach