

Solving Puzzles since 1946.

Creativity and Contribution



PHARMACEUTICAL NEWSLETTER

**Manufacture of Sulindac-Neusilin[®]
Amorphous Drug Complex Using Hot Melt Extrusion**

OVERVIEW

HOT MELT EXTRUSION (HME)

Hot Melt Extrusion (HME) offers several advantages over traditional pharmaceutical processing techniques such as fewer processing steps, continuous operation, little or no solvents, and commercial scale production of solid dispersions with improved bioavailability. New Chemical Entities (NCEs) with solubility problems are prime candidates for HME.

HME allows dispersion of poorly water-soluble drugs in a given matrix at the molecular level via the formation of a solid solution. The most exciting part of using **Neusilin®** in HME is that solid dispersions can be prepared without any addition of solvents or polymers.

In this newsletter, we report a summary from a recent publication by *Maclean et al.* on the *manufacture of Sulindac-Neusilin® amorphous drug complex via HME.*



Sulindac is a nonsteroidal anti-inflammatory drug (NSAID) used to treat mild to moderate pain and help relieve symptoms of arthritis.

HOT MELT EXTRUSION OF SULINDAC-NEUSILIN® DRUG COMPLEX


Blends of Sulindac-Neusilin® in 1:1 and 1:2 (w/w ratio) were prepared by mixing the components and charging them to the Brabender single-screw volumetric feeder that fed directly into the extruder hopper at a 5% feed rate.

The material was then extruded using a Prism PharmaLab 16 mm twin-screw extruder (25:1 L/D).

The screw speed was set to 50 rpm for the duration of the extrusion process. Conversion of the crystalline to amorphous complex was observed when HME was conducted at a temperature of 200°C, which is above the melting point of Sulindac.

The sample was recovered as powder from the HME apparatus. The dissolution profile of the amorphous powder complex is shown in Figure 1.

Figure 1: Dissolution profiles of amorphous complex and Crystalline Sulindac in 0.1N HCl




PHYSICAL / CHEMICAL STABILITY OF SULINDAC-NEUSILIN® HME COMPLEX

The HME samples remained amorphous after 3 months of storage at 40°C/75% RH (Figure 2). The samples were found to remain amorphous for more than one year at ambient conditions. The HPLC analysis during and after the storage period showed no chemical degradation thus confirming the physical and chemical stability of amorphous phase.

The stabilization of amorphous complex is an added advantage over other organic polymer excipients.



Figure 2: XRPD data of Sulindac-Neusilin® HME samples as a function of time at 40°C/75% RH.



All samples remained amorphous for 3 months at accelerated stress conditions.

This indicates that the Sulindac-Neusilin® HME samples maintained their physical stability and did not undergo crystallization or degradation during this time period.

SULINDAC-NEUSILIN® HME TABLETS

Approximately 100 g of 1:1 and 1:2 Sulindac-Neusilin® complexes were produced in this experiment.

Sulindac-Neusilin® complex powder recovered from HME was made into 100 mg, 200 mg and 500 mg tablets through direct compression.

Additional excipients and lubricants were added to match the commercial tablet while keeping the Sulindac weight constant (Table 1).

Tablets were made with crystalline Sulindac plus Neusilin® and HME complex to compare the dissolution profile.

Sulindac-Neusilin® 1:2 HME tablet showed 100% release in 90 minutes as against 9% release of Sulindac-Neusilin® crystalline tablets (Figure 3).






TABLE 1: COMPOSITIONS OF CRYSTALLINE AND HME AMORPHOUS 1:2 SULINDAC-NEUSILIN® TABLETS


Component	1:1 HME Tablet		1:2 Crystalline Tablet		1:2 HME Tablet	
	Weight (mg)	Weight (%)	Weight (mg)	Weight (%)	Weight (mg)	Weight (%)
Sulindac	200	40	50	25	50	25
Neusilin®	200	40	100	50	100	50
Starch 1500	47.5	9.5	20	10	20	10
Avicel PH 102	47.5	9.5	20	10	20	10
Stearic Acid	-	-	2	1	2	1
HPMC E4	-	-	8	4	8	4
Magnesium Stearate	5.0	1	-	-	-	-
Total Tablet Weight	500	100	200	100	200	100

Stabilizing mechanisms of ASD with Neusilin involve:

- Salt formation with acidic drugs like sulindac, indomethacin, ketoprofen, naproxen, ibuprofen and aceclofenac



- Hydrogen bonding



- Ion-dipole and dipole-dipole interactions for non-acidic drugs
- Moisture adsorption owing to high specific surface area thus preventing recrystallization of API





FIGURE 3: DISSOLUTION PROFILES OF HME SULINDAC-NEUSILIN® TABLETS



The tablets containing **Sulindac-Neusilin®** in a 1:2 ratio HME complex achieved complete drug release within 90 minutes. This indicates improved solubility and dissolution properties, potentially enhancing the drug's bioavailability and effectiveness.

CONCLUSION

HME was demonstrated as a method capable of large scale manufacturing of stable Sulindac- Neusilin® complex.

Sulindac-Neusilin® complex can be recovered as powder from HME apparatus.

No addition of organic polymers was necessary to prepare amorphous form of the drug.

HME samples remained amorphous after 3 months storage at 40°C/75% RH.

Both Physical and Chemical stability of the amorphous Sulindac-Neusilin® complex was maintained over 1 year at ambient temperature.

Tablets made from 1:2 Sulindac-Neusilin® amorphous complex showed 100 % release as against 9% release of crystalline drug.



Reference

Jenifer Maclean¹, Cesar Medina¹, Dominick Daurio¹, Fernando Alvarez-Nunez¹, Janan Jona¹, Eric Munson², Karthik Nagapudi¹

¹: Small Molecule Process and Product Development, Amgen, Thousand Oaks, California 91320

²: College of Pharmacy, University of Kentucky, Lexington, Kentucky 40536



**Fuji Chemical
Industries**

Solving Puzzles Since 1946.
Creativity and Contribution

CONTACT US

India

AstaReal (India) Private Limited
(A Fuji Chemical Group Company)
Unit No. 606/6th floor, Sunteck
Crest, Behind Mukund Hospital,
Near Leela Business Park,
Andheri East, Mumbai – 400059
Tel. +91-22-62369998
Email: pharma@fujichemical.co.jp

Europe

AstaReal AB
(A Fuji Chemical Group Company)
Forumvägen 14, Level 16,
131 53, Nacka, SWEDEN
Tel. +46-8-570-139-50
Email: pharma@fujichemical.co.jp

USA

Fuji Chemical Industries USA, Inc.
3 Terri Lane, Unit 12
Burlington, NJ 08016 USA
Tel. +1-609-386-3030
Email: contact@fujichemicalusa.com
www.fujichemicalusa.com

Japan

Fuji Chemical Industries Co., Ltd.
Shibakoen Ridge Building 2nd Floor,
1-8-21 Shiba Koen, Minato-ku,
Tokyo 105-0011 JAPAN
Tel. +81-3-3437-2350
Fax: +81-3-3437-2347
Email: pharma@fujichemical.co.jp
www.fujichemical.co.jp/english

*Regulations may vary from one region to another. Please contact
your local distributor or sales person for more information.