

PRODUCT SPECIFICATIONS SHEET

KLEPTOSE® HP PARENTERAL GRADE

Definition

Product Identifier

Product name: KLEPTOSE® HP PARENTERAL GRADE

HYDROXYPROPYLBETADEX is a partially substituted
poly (hydroxypropyl) ether of betadex.
INCI: HYDROXYPROPYL CYCLODEXTRIN

Specifications

A) CHARACTERS

| | |
|------------|---|
| APPEARANCE | White or almost white, amorphous or crystalline powder. |
| SOLUBILITY | Freely soluble in water and in propylene glycol. |

B) IDENTIFICATION

| | | |
|-----------------------|---------|----------|
| IDENTIFICATION-TEST A | EP / NF | COMPLIES |
| IDENTIFICATION-TEST B | EP / NF | COMPLIES |

C) TESTS

| | | |
|--|---------|---------------------|
| APPEARANCE IN SOLUTION | EP / NF | COMPLIES |
| CONDUCTIVITY | EP / NF | 200 micro.S/cm max. |
| RELATED SUBSTANCES | | |
| -IMPURITY A / BETADEX(on DS) | EP / NF | 1.5 % max. |
| -IMPURITY B / PROPYLENE GLYCOL(on DS) | EP / NF | 2.5 % max. |
| -SUM OF IMPURITIES OTHER THAN A(on DS) | EP | 1.0 % max. |
| -ANY OTHER IMPURITY(on DS) | NF | 0.25 % max. |
| -TOTAL OF OTHER IMPURITIES(on DS) | NF | 1 % max. |
| LOSS ON DRYING | EP / NF | 10.0 % max. |
| MOLAR SUBSTITUTION | EP / NF | 0.40 - 1.50 |
| PROPYLENE OXIDE(**) | NF | 0.0001 % max. |
| PARTICLE SIZE (Sieve) | | |
| - RESIDUE ON 315 microns | | 20 % max. |
| - RESIDUE ON 100 microns | | 50 % min. |
| MICROBIOLOGICAL VALUES: | | |
| -TOTAL AEROBIC MICROBIAL COUNT | EP / NF | 100 CFU/g max. |

PRODUCT SPECIFICATIONS SHEET

KLEPTOSE® HP PARENTERAL GRADE

| | | |
|--------------------------------|---------|---------------------|
| -TOTAL YEASTS AND MOULDS COUNT | EP / NF | 100 CFU/g max. |
| -ESCHERICHIA COLI | EP | Not detected in 10g |
| -SALMONELLA | EP | Not detected in 10g |
| BIOLOGICAL VALUES: | | |
| - BACTERIAL ENDOTOXINS | EP | <10 EU/g |

Indicatives Values

| | | |
|----------------------------------|---------|-------------|
| MOLAR SUBSTITUTION NOMINAL VALUE | EP / NF | 0.90 |
| MOLAR SUBSTITUTION (MS) | EP / NF | 0.81 - 0.99 |

Comments

This substance is suitable for use in the manufacture of parenteral preparations.

Methods used by Roquette may be the Pharmacopoeia methods or alternative validated methods which have been compared to the Pharmacopoeia methods.

Caption

- "EP" stands for European Pharmacopoeia
- "NF" stands for National Formulary from USP-NF
- (**) Monitoring plan

Conformity

Conforms to the requirements of the current monograph

- **European Pharmacopoeia (EP)** HYDROXYPROPYLBETADEX (1804)
- **National Formulary from USP-NF** HYDROXYPROPYL BETADEX

Please contact us for any statement regarding compliance to the General Chapters (elemental impurities, residual solvents, organic volatile impurities, metal catalyst, metal reagent).

Storage

RETEST DATE Manufacturing date + 3 years, in its unopened packaging.

- The product durability may vary according to packaging type and manufacturing plant. Proper information is shown on labelling and CoA.
- We recommend to preserve the product in its unopened original packaging, preferably protected from wide variations of temperature and humidity.
- Upon opening, use the product as quickly as possible to prevent moisture regain.

KLEPTOSE® HP PARENTERAL GRADE

Disclaimer

The information provided in this Product Specification Sheet relates only to the specific product designated and may not be applicable when such product is used in combination with other materials or in any process.

All information and instructions provided in this Product Specification Sheet are based on the current state of our knowledge at the latest revision date indicated. It is the responsibility of the user to be aware of and to follow the regulations applying to our product for its possession, handling and use.

Notes : All the dates are formatted like YYYY/MM/DD.